

A 3D anatomical illustration of a human head and neck in profile, facing left. The brain and spinal cord are rendered in a vibrant blue color, while the rest of the head and neck are shown in a lighter, semi-transparent blue. The background is a light blue gradient with a pattern of overlapping circles of varying sizes.

Advanced Neuromonitoring Solutions

**Multimodal Neuromonitoring
with Measuring Catheters**

Table of Contents

- 3** Multi-Modal Neuromonitoring
- 5** Extensive Measuring Catheter Portfolio
- 6** NEUROVENT Benefits
- 8** Transferring Measurement Values to the Patient Monitor
- 10** Continuous Measurement of Intracranial Pressure
- 12** NEUROVENT-PTO
- 14** Touch-screen Monitors
- 18** Data Display Devices
- 20** Software Solution
- 22** Application of Catheters
- 24** Product List
- 31** What can we do for you?



Measuring Catheters

Important parameters measured with high-precision microchip catheters

- **ICP (intracranial pressure)**
- **ICT (intracranial temperature)**
- **p_{ti}O₂ (oxygen partial pressure)**

The multi-modal neuromonitoring performed with the measuring catheter in the field of neurosurgery enables early recognition of potential cerebral damages.

ICP is measured using semiconductor pressure sensors. The quenching process of fluorescence is used to measure p_{ti}O₂. Consequently, the level and changes in the parameters are measured safely, quickly and accurately.





- Epidural
- Parenchymal
- Ventricular

Measuring Locations

Parenchymal

- **NEUROVENT-P / NEUROVENT-PX**
Parenchymal ICP measurement
- **NEUROVENT-P-TEMP**
Parenchymal ICP and temperature measurement
- **NEUROVENT-PTO Series**
Parenchymal ICP, temperature and p_vO₂ measurement
- **NEUROVENT-TO**
Parenchymal temperature and p_vO₂ measurement



Ventricular

- **NEUROVENT Series**
Ventricular ICP measurement and CSF-Drainage
- **NEUROVENT-TEMP Series**
Ventricular ICP and temperature measurement with CSF-Drainage
- **NEUROVENT-Sleeve Housing**
Ventricular CSF-Drainage and parenchymal ICP measurement
- **NEUROVENT VP 16**
Ventricular CSF-Drainage and parenchymal ICP measurement, neuronavigable



Epidural

- **NEURODUR**
Epidural ICP measurement
- **NEURODUR-TEMP**
Epidural ICP and temperature measurement



Clinical Advantages


of RAUMEDIC NEUROVENT Catheters

- + Plug & Play system – no catheter calibration required
- + Compatible with all standard patient monitors
- + Excellent measurement stability and linearity
- + MR conditional at 1,5 T and 3,0 T¹ – no surgical intervention and disposition of the catheter required
- + Monitor change without measurement loss of ICP is possible using zero point simulator (NPS2)*

¹ Indicated within non-clinical laboratory tests.
*For more information about the NPS2 please have a look at page 25.

Comparison of Material Surfaces

RAUMEDIC ventricular catheters are made from polyurethane (PU) – compared to standard ventricular catheters that are typically made from silicone. The difference of the distinct catheters' surface structures becomes clear under a scanning electron microscope (SEM).



» Innovative Materials » Smooth Surface Structure



Advantages of RAUMEDIC microchip catheter NEUROVENT in comparison with an EVD system

- Continuously ICP measurement and CSF drainage
- Precise measurement of the pressure variation
- No artefacts by moving the patient or opening of the CSF-Drainage valve
- No hydrostatic failure influences
- Fast adaption of changes in ICP



8 Accessories

The Zero-Point Simulator NPS2 is attached directly to the invasive blood pressure (IBP) port of the patient monitor.

All RAUMEDIC catheters are calibrated in the manufacturing process, and therefore are ready to use.

The Zero-Point Simulator NPS2 is used to transfer the ,zero' to the patient monitor.


Your advantages

- Adapters available for all common patient monitors*
- Direct connection to the patient monitor
- No ICP monitor required
- Plug & Play system – no catheter calibration required



*Please find the full list of adapters as well as further technical data and product information on page 25.

9 Measuring chain from the catheter to the patient monitor



Please find further technical data and product information on page 25.



Measuring ICP during Patient Transport



Measuring chain from the catheter to the NPS3 pressure display



The intrahospital transportation of critically ill patients with severe brain diseases is linked to a considerable rate of complications. Therefore, performing the transport of critically ill under monitoring ICP can often be crucial for positive patient outcomes.


The battery-powered pressure display device NPS3 can simply be connected to our ICP-TEMP-Cable, which – in turn – can be connected to a measuring catheter.

Great Advantages during Patient Transport

- Monitoring of ICP during patient transport
- Battery-powered display device
- No further monitor required
- Plug & Play system – no catheter calibration required



One Catheter, three Measurement Functions



For measuring $p_{ti}O_2$, the quenching process of fluorescence is used. At the same time, parenchymal pressure is measured with the same measuring catheter via semiconductor pressure sensors.

Oxygen partial pressure measurement records the available oxygen in the brain tissue. This ensures possible cerebral damage is quickly detected and appropriate measures for the avoidance of cerebral ischaemia can be taken.

Unique advantages of a unique product



- Unique catheter that measures three parameters at once
- Parenchymal pressure, temperature and $p_{ti}O_2$ measurement in one catheter
- Easy handling via Plug & Play system – no calibration required
- No oxygen consumption by the O_2 sensor
- No refrigeration required
- Data recording and data display devices available*



NEUROVENT-PTO


NEUROVENT-TO

Catheter for measuring temperature and $p_{ti}O_2$



NEUROVENT-PTO 2L

Specially developed catheter for craniotomies which measures ICP, temperature and $p_{ti}O_2$



NEUROVENT-PTO 2L BOLT

Catheter for measuring ICP, temperature and $p_{ti}O_2$ for application with BOLT KIT PTO 2L and a microdialysis catheter



BOLT KIT PTO 2L

Two lumen BOLT for safe and functional implantation of the NEUROVENT-PTO 2L BOLT and a microdialysis catheter



*Please find further technical data and product information on page 26.



Smart Neuro-monitoring


The Next Generation of Smart Neuromonitoring Devices

For the calculation and visualization of vital parameters

The RAUMED NeuroSmart is available for the visualization and storage of ICP and ICPT (telemetrically) measurement data. In another version, the RAUMED NeuroSmart logO, the oxygen partial pressure $p_{t_i}O_2$ can also be recorded and visualized.



RAUMED NeuroSmart logO



Display of ORx, PRx and CPPopt
Regression curve CPPopt

Live data
Live Graph View

Advantages and Features

Features

- Colored touch-screen display
- Colors can be set individually
- Calculation of PRx, ORx and CPPopt
- Various scaling options
- Display of the curves and trend graph
- 2 USB interfaces (USB stick and PC connection)

Clinical Advantages

- Display of ICP, brain temperature, $p_{t_i}O_2$, ART, CVP, ICPA-amplitudes
- Integrated data storage for up to 10 days
- Audible and visual alarms
- Battery / mains operation possible
- Mobile use
- Connection to the patient monitor possible
- Attachment to the pole



Device features

- Invasive pressure (2 x)
- Telemetry pressure (1 x)
- Oxygen partial pressure (1 x) (only RAUMED NeuroSmart logO)
- Temperature (2 x ICT)
- Analogue outputs (2 x)
- USB interfaces (2 x)
- Analogue Rec output (1 x)

Connections RAUMED NeuroSmart logO



2 USB interfaces

General information

| | |
|---|--|
| Display | LCD, color, 10-inch |
| Graphs / Trend display | Selectable via menu |
| Alarm limits | Selectable via menu |
| Dimensions | Approx. 310 x 225 x 150 mm (W x L x D) |
| Mass | Approx. 3 kg with battery and stand holder |
| Power supply | Internal battery with power adapter |
| Operating period in battery mode | ≥ 4 h (charged, new battery, all channels occupied) |



More Monitoring Solutions

Keeping an eye on key vital functions at all times

In neurocritical or neurointensive care, the monitoring of vital physiological functions of patients plays a central role. The MPR2 logO DATALOGGER and EASY logO from RAUMEDIC are comprehensive monitoring solutions for this.


EASY logO: Simple display of data

For displaying of ICP data, temperature, and oxygen partial pressure.



MPR2 logO DATALOGGER: Easy visualizing and recording

For displaying and recording ICP, temperature and oxygen partial pressure – and visualization of parameters as curves and trend graphs.




| | | |
|--|---|---|
| Mains operation | ✓ | ✓ |
| Rechargeable battery | ✗ | ✓ |
| 2 x Analog outputs (Transfer of pressure value) | ✓ | ✓ |
| USB interface | ✗ | ✓ |
| Data storage | ✗ | ✓ |
| Curve display | ✗ | ✓ |
| Possible display of | | |
| • ICP | ✓ | ✓ |
| • p _t O ₂ | ✓ | ✓ |
| • Temperature | ✓ | ✓ |
| • ICPA | ✓ | ✓ |
| • ART | ✗ | ✓ |
| • CPP | ✗ | ✓ |



RAUMED DataView

Software for transferring measurement data to a PC or Laptop

RAUMED DataView is a software* solution for transferring recorded and saved data on the RAUMEDIC devices to a computer for visualization.



*No medical product according to Regulation (EU) 2017/745. The application is not required for the operation and intended use of the devices.

»Benefits of RAUMED DataView

- + No internet connection required
- + Display of ICP, temperature, $p_{ti}O_2$ and ART on the PC
- + Extended display of PRx, ORx and CPPopt
- + Graphic (comparative) representation of data
- + Filter functions for data search
- + Data export to CSV, EDF, PDF and RAUMED DataView format (dv.data)



Application Accessories


Single-use application accessories for the subcutaneous tunneling of RAUMEDIC catheters.

Spliceable Tunneling Sleeve

RAUMEDIC Tunneling Sleeves are made of biocompatible, polymer material – in-vitro tested according to EN ISO 10993-1.

The sleeve is chamfered for low resistance application.

Available in two versions: for parenchymal and ventricular catheters, NEUROVENT-PTO 2L.



Tunneling KIT


With its fir-tree-like connection of the trocar, the Tunneling KIT provides a secure hold for the tunneling sleeve, which is cut off after being placed under the scalp.

Available in two versions: for parenchymal and ventricular catheters, NEUROVENT-PTO 2L.



BOLT-DRILL KIT

A recognized method for inserting a catheter is to use a drill and an appropriate drill bit to pass it through the calvaria. To safely position the measuring catheter, a BOLT is screwed into the drilled hole. Our BOLT-DRILL KIT offers the necessary components in a single set.



BOLT KIT components

- 1. Polymer screw with fixing cap and sealing ring
- 2. Dura opener
- 3. Screw-in tool

DRILL KIT components

- 4. Drill bit with stopper
- 5. Allen key

BOLT KIT advantages

- Material suitable for all imaging methods
- Low BOLT height
- Self-cutting thread with sealing function



Product List

Parenchymal measurements

| Product | Version | Dimension | Article number |
|-----------------------|---|-----------|----------------|
| NEUROVENT-P | ICP | 5F | 092 946-001 |
| NEUROVENT-PX | ICP | 5F | 091 580-001 |
| NEUROVENT-P-TEMP | ICP + temperature | 5F | 094 268-001 |
| NEUROVENT-PTO | ICP + temperature + p _i O ₂ | 5F | 095 008-001 |
| | application with BOLT-DRILL KIT PTO | | 092 380-001 |
| NEUROVENT-PTO 2L | ICP + temperature + p _i O ₂ | 5F | 095 108-001 |
| NEUROVENT-PTO 2L BOLT | ICP + temperature + p _i O ₂ | 5F | 095 308-001 |
| | application with BOLT KIT PTO 2L | | 096 076-001 |
| | and DRILL KIT CH9 | | 091 668-002 |
| NEUROVENT-TO | Temperature + p _i O ₂ | 3F | 095 908-001 |
| | application with BOLT-DRILL KIT PTO | | 092 380-001 |

Ventricular measurements

| Product | Version | Dimension | Article number |
|--------------------------|---|-----------|----------------|
| NEUROVENT | ICP + drainage, with stylet | 9F | 092 956-001 |
| NEUROVENT 6F | ICP + drainage | 6F | 094 678-001 |
| NEUROVENT-IFD-S | ICP + drainage, soft internal guide wire | 9F | 091 678-001 |
| NEUROVENT-IFD-R | ICP + drainage, rigid internal guide wire | 9F | 095 317-001 |
| NEUROVENT-Sleeve Housing | ICP (parenchyma), drainage (ventricular), with sleeve housing | 9F | 091 576-001 |
| NEUROVENT VP 16 | ICP (parenchyma), drainage (ventricular), neuronavigable | 9F | 096 704-001 |
| NEUROVENT-TEMP | ICP + drainage + temperature, with stylet | 9F | 094 278-001 |
| NEUROVENT-TEMP-IFD-S | ICP + drainage + temperature soft internal guide wire | 9F | 094 288-001 |
| NEUROVENT-TEMP-IFD-R | ICP + drainage + temperature rigid internal guide wire | 9F | 095 327-001 |

Epidural measurements

| Product | Version | Dimension | Article number |
|---------------|-------------------|------------------------------------|----------------|
| NEURODUR | ICP | 5.8 x 2.1 mm (measurement head) | 092 976-001 |
| NEURODUR-TEMP | ICP + temperature | 5.8 x 2.1 mm (measurement head) | 094 298-001 |

Catheters

Technical data

| | |
|---|--------------------------------------|
| Pressure measurement range | -40 to +400 mmHg (-5.3 to 53 kPa) |
| Upper cut-off frequency | 20,000 Hz (-3 dB) |
| Catheter material | Polyurethane |
| Measurement range temperature sensor | +25°C to +45°C |
| Pressure sensitivity | 5 µV/V/mmHg |
| Measurement range p _i O ₂ * | 0-200 mmHg |
| Electrical catheter length (tip to connector) | |
| - Parenchymal | approx. 55 cm |
| - Ventricular | approx. 55 cm |
| - Epidural | approx. 55 cm |

* Measurement accuracy ± 2.5 mmHg (for < 120 mmHg)

Zero Drift Pressure

Ø Deviation 0.6 mmHg after 5 days*

* Bench test assessment of the new Raumedic Neurovent-P ICP sensor: a technical report by the BrainIT group Citerio G., Piper I., Cormio M., Galli D., Cazzaniga S., Enblad P., Nilsson P., Contant C., and Chambers I., BrainIT Group Acta Neurochirurgica (Wien). 2004, Aug; DOI: 10.1007/s00701-004-0351-z



Connecting Cables

| Product | Description | Article number |
|------------------------------|--|----------------|
| ICP-TEMP-Cable | Connecting cable between ICP catheter and zero point simulator NPS2 | 094 328-001 |
| ICP-TEMP-Adapter | Adapter between zero point simulator NPS2 and patient monitor | 094 323-001 |
| ICP-TEMP-Adapter Philips/HP | Adapter between zero point simulator NPS2 and patient monitor Philips/HP | 094 047-001 |
| NPS2 Siemens/Dräger Infinity | Adapter cable to Siemens/Dräger Infinity patient monitor | 092 627-001 |
| NPS2 Philips/HP | Adapter cable to Philips/HP patient monitor | 092 637-001 |
| NPS2 Nihon Kohden BSM 41xx | Adapter cable to Nihon Kohden BSM 41xx patient monitor | 094 716-001 |
| NPS2 GE/MARQUETTE | Adapter cable to GE/MARQUETTE patient monitor | 093 807-001 |
| NPS2 GE | Adapter cable to GE patient monitor | 093 999-001 |
| NPS2 SpaceLabs | Adapter cable to SpaceLabs patient monitor | 091 715-001 |
| NPS2 Fukuda Denshi | Adapter cable to Fukuda Denshi patient monitor | 096 003-001 |
| NPS3 | Battery operated pressure device | 091 656-001 |

Zero point simulator NPS2 for further monitor types upon request

Oxygen partial pressure measurement

| Product | Version | Dimension | Article number |
|------------------------|--|-----------|----------------|
| NEUROVENT-PTO | ICP + temperature + p _t O ₂ application with BOLT-DRILL KIT PTO | 5F | 095 008-001 |
| | | | 092 380-001 |
| NEUROVENT-TO | Temperature + p _t O ₂ application with BOLT-DRILL KIT PTO | 3F | 095 908-001 |
| | | | 092 380-001 |
| NEUROVENT-PTO 2L | ICP + temperature + p _t O ₂ | 5F | 095 108-001 |
| NEUROVENT-PTO 2L BOLT | ICP + temperature + p _t O ₂ application with BOLT KIT PTO 2L and DRILL KIT CH9 | 5F | 095 308-001 |
| | | | 096 076-001 |
| | | | 091 668-002 |
| EASY logO | Data display | | 095 264-002 |
| MPR2 logO DATALOGGER | Data recording and storage | | 095 254-002 |
| RAUMED NeuroSmart logO | Data recording and storage | | 095 294-001 |

RAUMED NeuroSmart and Accessories

| Product | Article number |
|---------------------------|----------------|
| RAUMED NeuroSmart | 095 284-001 |
| ICP-TEMP-Cable | 094 328-001 |
| RAUMED DataView* | 296 900-001 |
| USB-Cable* | 283 949-001 |
| Power adapter NeuroSmart | 284 037-001 |
| Rec-BNC-Cable NeuroSmart* | 096 096-001 |



RAUMED NeuroSmart logO and Accessories

| Product | Article number |
|---------------------------|----------------|
| RAUMED NeuroSmart logO | 095 294-001 |
| Cable PTO | 095 624-001 |
| Cable LWL | 095 657-001 |
| ICP-TEMP-Cable | 094 328-001 |
| RAUMED DataView* | 296 900-001 |
| USB-Cable* | 283 949-001 |
| Power adapter NeuroSmart | 284 037-001 |
| Rec-BNC-Cable NeuroSmart* | 096 096-001 |

*No medical product according to Regulation (EU) 2017/745.



Accessories for all devices

| Product | Article number | RAUMED NeuroSmart | RAUMED NeuroSmart logO | MPR2 logO DATALOGGER | EASY logO |
|----------------------------------|----------------|-------------------|------------------------|----------------------|-----------|
| Cable PTO | 095 624-001 | x | x | x | |
| Cable LWL | 095 657-001 | | x | x | x |
| ICP-TEMP-Cable | 094 328-001 | x | x | x | |
| ICP-TEMP-Adapter | 094 323-001 | | | x | |
| Main power adapter EASY logO | 284 017-001 | | | | x |
| Wide range power adapter MPR 1/2 | 284 027-001 | | | x | |
| Stand holder DATALOGGER | 283 957-002 | | | x | x |
| Table stand DATALOGGER | 283 959-002 | | | x | x |
| RAUMED DataView* | 296 900-001 | x | x | x | |
| USB-Cable* | 283 949-001 | x | x | x | |
| Power adapter NeuroSmart | 284 037-001 | x | x | | |
| Rec-BNC-Cable NeuroSmart* | 096 096-001 | x | x | | |



Connecting cables from RAUMEDIC device to patient monitor

| Product | Article number | GE/MARQUETTE | Philips/HP | Siemens/Dräger Infinity | SpaceLabs | Nihon Kohden 41xx |
|------------------|----------------|--------------|------------|-------------------------|-----------|-------------------|
| Cable DATALOGGER | 094 858-001 | x | x | x | x | |
| Cable DATALOGGER | 094 868-002 | x | x | x | x | |
| Cable DATALOGGER | 094 878-002 | x | x | x | x | |
| Cable DATALOGGER | 094 967-001 | x | x | x | x | |
| Cable DATALOGGER | 095 017-001 | x | x | x | x | |



Transducercables between RAUMEDIC device and disposable transducer

| Product | Article number | Medex MX 960 | Edwards TRUWAVE | Becton Dickinson | Combitrans | pvb xtrans |
|-----------------|----------------|--------------|-----------------|------------------|------------|------------|
| Transducercable | 095 974-001 | x | x | x | | |
| Transducercable | 096 036-001 | x | x | x | | |
| Transducercable | 096 046-001 | x | x | x | | |
| Transducercable | 096 664-001 | x | x | x | | |
| Transducercable | 096 494-001 | x | x | x | | |



Accessories for Catheters

Spliceable Tunneling Sleeve

| Product | Article number |
|---|----------------|
| Spliceable Tunneling Sleeve CH8 (for parenchymal catheters) | 090 506-002 |
| Spliceable Tunneling Sleeve CH12 (for ventricular catheters and NEUROVENT-PTO 2L) | 090 717-001 |

Tunneling KIT

| Product | Article number |
|---|----------------|
| Tunneling KIT CH8 (for parenchymal catheters) | 090 516-001 |
| Tunneling KIT CH12 (for ventricular catheters and NEUROVENT-PTO 2L) | 090 727-001 |



BOLT-DRILL KIT

| Product | Version | Article number |
|----------------------|--|----------------|
| BOLT KIT CH5 | For parenchymal catheters | 091 868-002 |
| DRILL KIT CH5 | For BOLT KIT CH5 | 091 878-002 |
| BOLT-DRILL KIT CH5 | Set for parenchymal catheters | 091 888-001 |
| BOLT KIT CH9 | For ventricular catheters | 091 688-002 |
| DRILL KIT CH9 | For BOLT KIT CH9 | 091 668-002 |
| BOLT-DRILL KIT CH9 | Set for ventricular catheters | 091 898-001 |
| BOLT KIT PTO | Only for NEUROVENT-PTO/-TO | 096 026-001 |
| BOLT-DRILL KIT PTO | Set for NEUROVENT-PTO/-TO | 092 380-001 |
| BOLT-DRILL KIT VP 16 | Only for NEUROVENT VP 16 and NEUROVENT-Sleeve Housing | 092 969-001 |
| RALK-Hand Drill | Autoclavable drill | 231 584-002 |



30 References

Poster (2015) Medstar Washington Hospital Center, Washington, D.C., RAUMEDIC Bolt: Initial clinical experience in a neurosurgical population, MD Rocco Armonda, MD Daniel Felbaum, MD Kyle Mueller, MD Anthony Conte, MD R. Bryan Mason, MD Edward Aulisi;

Journal of Clinical Neuroscience (2011), DOI:10.1016/j.jocn.2011.04.026, An outcome analysis of two different procedures of burr-hole trephine and external ventricular drainage in acute hydrocephalus, Petra Schödel, Martin Proescholdt, Odo-Winfried Ullrich, Alexander Brawanski, Karl-Michael Schebesch;

www.neurosurgery-online.com (2010), Neurosurgery 67:1716-1723, Evaluation of a Novel Brain Tissue Oxygenation Probe in an Experimental Swine Model, MD Berk Orakcioglu, MD Oliver W. Sakowitz, MD Jan-Oliver Neumann, MD Modar M. Kentar, MD PhD Andreas Unterberg, MD PhD Karl L. Kiening;

Acta Neurochir (2009) DOI 10.1007/s00701-009-0532-x, Brain tissue oxygen monitoring: a study of in vitro accuracy and stability of NEUROVENT-PTO and Licox sensors, Karlis Purins, Per Enblad, Bo Sandhagen, Anders Lewén;

Acta Neurochir (Wien) (2004) DOI 10.1007/s00701-004-0351-z, Bench test assessment of the new RAUMEDIC NEUROVENT-P ICP sensor: a technical report by the BrainIT group, G. Citerio, I. Piper, M. Cormio, D. Galli, S. Cazzaniga, P. Enblad, P. Nilsson, C. Contant, and I. Chambers on behalf of the BrainIT Group;

Journal of Neuroscience Methods 139 (2004) 161-165, Accuracy and stability of temperature probes for intracranial application, Beat Alessandri, Bernd M. Hoelper, Robert Behr, Oliver Kempfski;

Acta Neurochir (2003) 145: 185-193, DOI 10.1007/s00701-002-1052-0, Clinical evaluation of a new intracranial pressure monitoring device, R. Stendel, J. Heidenreich, A. Schilling, R. Akhavan-Sigari, R. Kurth, T. Picht, T. Pietilä, O. Suess, C. Kern, J. Meisel, and M. Brock.

What can we do for you?

Global Service

31

Location | RAUMEDIC Group

German Headquarters

RAUMEDIC AG
Hermann-Staudinger-Str. 2
95233 Helmbrechts
T +49 9252 359-1587
F +49 9252 359-513333 (orders)

neuromonitoring@raumedic.com
raumedic.com/neuromonitoring

